Michelson Experiment Summary

Yongao Hu

MIT Department of Physics

February 26, 2025

Outline

Introduction

- Light exhibits interference when two waves overlap.
- The Michelson interferometer, designed by Albert Michelson, measures the wavelength of light.
- Historically, it disproved the existence of the "luminiferous aether" and contributed to relativity.
- Interferometers are used in:
 - Precision measurements
 - LIGO for detecting gravitational waves [?]



Interference of Light Waves

• The interference pattern depends on the phase difference:

$$\Delta \phi = 2\pi \frac{\Delta x}{\lambda}$$

• Intensity of interference:

$$I = I_0 \cos^2(\Delta \phi)$$

- Maximum interference: $\Delta \phi = 2\pi n$
- Minimum interference: $\Delta \phi = (2n+1)\pi$



Visibility of the Interference Pattern

• Visibility measures how well the pattern is resolved:

$$V = \frac{I_{\text{max}} - I_{\text{min}}}{I_{\text{max}} + I_{\text{min}}}$$

- Ideal case: V=1
- ullet Low visibility \Rightarrow poor alignment or coherence issues.



Piezoelectric Transducer for Path Difference

- A piezoelectric transducer moves a mirror in response to voltage.
- This changes the path length Δx , altering the interference pattern.
- The wavelength of light is found using:


$$\lambda = 4\Delta V \times \frac{\Delta L}{V}$$



Michelson Interferometer Setup

- Light source: laser
- \bullet 50/50 beam splitter splits the light into two paths.
- One mirror is stationary, the other mounted on a piezoelectric transducer.
- Detector: photodiode (voltage output proportional to intensity).
- Signal generator applies a triangular voltage to move the mirror.


Experimental Setup - Schematic



Data Processing

- Voltage signals recorded from signal generator and photodiode.
- Used findpeak in SciPy to locate primary maxima and minima.
- Applied KMeans clustering (SciKit-Learn) to group voltage values.
- Measured ΔV between consecutive maxima and minima.

Interference Pattern Visibility



Visibility Measurement

- ullet Maximum voltage: (1.592 ± 0.002) V
- ullet Minimum voltage: (1.110 ± 0.002) V
- Visibility:

$$V = 0.179 \pm 0.001$$

• Imperfect alignment likely lowered visibility.



Wavelength Measurement


- ullet Measured voltage difference between peaks: (3.192 \pm 0.055) V
- Calculated wavelength:

$$\lambda = 601 \pm 37 (ext{sys}) \pm 10 (ext{stat}) \; ext{nm}$$

Matches expected wavelength of the orange laser.



Voltage Data at Peaks



Sources of Uncertainty

• **Systematic error**: Calibration of the piezoelectric device

$$\delta\lambda=$$
 37 nm

• **Statistical error**: Peak detection noise

$$\delta\lambda=$$
 10 nm

• **Total error dominated by systematic effects.**

Conclusion

• Successfully measured the wavelength of an orange laser:

$$601 \pm 37 (ext{sys}) \pm 10 (ext{stat})$$
 nm

- Demonstrated interference and wave nature of light.
- Future improvements:
 - Improve mirror alignment for higher visibility.
 - Use a more precise piezoelectric transducer.
 - Implement automated feedback to stabilize fringes.



Acknowledgments

- Special thanks to J. Lewis, M. Bohdan, and V. Tran for assistance.
- Thanks to the MIT 8.13 teaching team for guidance.
- Supported by the MIT Department of Physics.

References I

